

Quick-RNATM MiniPrep

- Sample Lysis
 - Remove media and lyse cells in 300 μ L (up to 5x10⁶ cells) or 600 μ L (>5x10⁶ cells) RNA Lysis Buffer by scraping
 - Put in 1.5mL tub
- Sample Clearing and gDNA Removal
 - Centrifuge lysate at 13,000g for 1min
 - Transfer the supernatant into a **Spin-Away Filter** in a Collection Tube and centrifuge at 13,000g for 1min
 - Save the flow-through
- RNA Purification
 - Add 1 volume ethanol (95-100%) to the sample in **RNA Lysis Buffer** (1:1). Mix well.
 - Transfer the mixture to a **Zymo-Spin IICG Column** in a Collection Tube and centrifuge for 30 seconds. Discard the flow-through
 - **In-column DNase I Treatment** (optional)
 - Add 400 μ L **RNA Wash Buffer** to the column and centrifuge for 30sec. Discard the flow-through
 - Add 80 μ L **DNase I reaction mix** directly to the column matrix. Incubate at room temp for 15min, then centrifuge for 30 sec

For each sample:

DNase I	5μL
10X DNase I Reaction Bffr	8μL
DNase/RNase-Free Water	3μL
RNA Wash Bffr	64μL

 - Add 400 μ L **RNA Prep Buffer** to the column and centrifuge for 30 sec, discard the flow-through
 - Add 700 μ L **RNA Wash Buffer** to the column and centrifuge for 30 sec, discard the flow-through
 - Add 400 μ L **RNA Wash Buffer** and centrifuge the column for 2min
 - Place the column in RNase-free tube. Add 30 μ L **RNase-Free Water** directly to the column, centrifuge at 16000g for 30sec
 - Nanodrop RNA