

## Packaging of shRNA constructs into Viral Particles and Infection of Target Cells

### Day 1

1. Plate 293FT cells in 4mL DMEM/10%FBS on a 60mm plate ( $1 \times 10^6$  cells). Prepare one 60mm plate for each construct that must be packaged into viral particles. Aliquots of 293FT cells are stored in liquid nitrogen Rack 3 Box 11

### Day 2

2. Begin transfection ~24h later. For each transfection, place 200 $\mu$ l sterile 150mM NaCl into a 1.5mL eppendorf tube.
3. Add 0.5 $\mu$ g of each packaging vector (PI, PII, PIII or pCMV-VSV-G, p $\Delta$ 8.9 ). Do not premix packaging vectors.
4. Add 1 $\mu$ g lentiviral plasmid and gently mix
5. Add 16 $\mu$ l TurboFect. Gently mix and incubate at RT for 15 min.
6. Add to 293FT cells by gently dripping into the media. Incubate overnight at 37°C.

### Day 3

7. In the morning, aspirate media from the 293 FT cells and replace with 4 mL of appropriate target media. Be careful not to disturb the 293FT cells as they are loosely adherent.
8. Plate target cells for infection tomorrow.

### Day 4

9. Prior to infection; pretreat target cells with polybrene (Sigma H9268, stored in Media Fridge). 1000x stock prepared in sterile water is stored at 4°C. Add 1 $\mu$ l of stock per mL of target cell media (10 $\mu$ l for 10 mL media). Return to incubator for 30-60 minutes.
10. Pipet viral media from the 293FT cells into 15mL conical tubes. Add back 4mL media to 293FT packaging cells.
11. Add 1 $\mu$ l polybrene stock per mL viral media (total 4 $\mu$ l for 4mL viral media). Filter through 0.45 $\mu$ m syringe filter.
12. Aspirate media from target cells
13. Split viral media equally among target cells and return to incubator for 1 hr. For cells in suspension, resuspend  $1 \times 10^7$  cells in 200 $\mu$ l viral media and spin at 2000rpm for 1hr at RT.
14. Aspirate viral media (off of target cells) and replace with normal target cell media (10mL). Incubate at 37°C overnight. For cells that are difficult to infect, add an additional volume of polybrene-containing media and incubate for an additional 1hr or up to overnight

### Day 5

15. Repeat day 4. (This greatly increases transfection efficiency)

### Day 6

16. Aspirate target cell media and replace with fresh media containing 2 $\mu$ g/mL puromycin to begin positive selection of infected cells. Puromycin stock is 1000x and stored at -20°C.

### Days 9-13

17. Observe puromycin resistant colonies. Feed and expand as necessary.
18. Screen for knockdown by western blot and PCR. Begin cloning of successful knockdowns
19. Freeze aliquots for storage